

USP 467 Residual Solvents (OVIs) Headspace on MEGA-624 FAST (G43 Phase)

Column: **MEGA-624 FAST - 0.10mm, 0.45µm, 10m**

Catalog Code: F-624-010-045-10

Conditions:

Injection: Split 250°C, 0.5mL with Gas Syringe, 1:100 Split Ratio.

Sample: Headspace of a Residual Solvent Mix hold 45min @ 80°C.


Detector: FID 250°C.

Oven Temp.: 35°C, 15°C/min, 100°C.

Carrier Gas: Hydrogen, 0.4 mL/min.

Peaks:

1. Methylene Chloride
2. Chloroform
3. Benzene
4. 1,1,1-trichloroethylene
5. 1,4-dioxane



Analysis carried out with DANI MASTER-GC.

