



BIOLOGIX社 Filter Chip Compatibility sheet

Pipette Brand		Eppendorf Research Plus						Gilson Pipetman(PG)						Gilson Pipetman(PN)						Gilson Pipetman(PL)					
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	10-100ul	20-200ul	100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul
Volume	Product No.	Universal Fit Filter Tips																							
0.5-10ul	23-0010S	√	√					√	√					√	√					√	√				
	23-0011S	√	√					√	√					√	√					√	√				
1-20ul	23-0020S			√						√												√			
	23-0021S			√						√												√			
1-50ul	23-0050S				√						√	√					√	√					√	√	
	23-0051S				√						√	√					√	√					√	√	
1-100ul	23-0100S				√						√	√					√	√					√	√	
	23-0101S				√						√	√					√	√					√	√	
1-300ul	23-0200S					√						√						√						√	
	23-0201S					√						√						√						√	
100-1000ul	23-1000S						√					√						√						√	
	23-1001S						√					√						√						√	

