

Amines - Basic Compounds

Column: **MEGA-PS255 - 0.25mm, 3.00µm, 30m**

Catalog Code: C-PS255-025-300-30

Conditions:

Injection: Split 250°C.

Detector: FID 300°C.

Oven Temp.: 100°C, 2.5°C/min, 200°C.

Carrier Gas: Hydrogen, 80kPa.

Peaks:

1. 2-Dimethylamino ethanol
2. 1,4-Dimethyl piperazine
3. Diethanolamine
4. Methyl-diethanolamine

