



## Organic Acids

Column: **MEGA-ACID (FFAP) - 0.53mm, 1.00 $\mu$ m, 15m**

Catalog Code: C-ACID-053-100-15

### Conditions:

Injection: Split 250°C, 0.2 $\mu$ L, 1:140 Split Ratio.

Detector: FID 250°C.

Oven Temp.: 200°C Isothermal.

Carrier Gas: Helium, 20 kPa.

### Peaks:

1. Sorbic acid
2. Dehydroacetic acid
3. Benzoic acid

