

Peaks:

- | | | |
|----------------|------------------------|---------------------|
| 1. Trichlorfon | 8. Chlorpyrifos methyl | 15. Quinalphos |
| 2. Phosdrin | 9. Parathion | 16. Bromophos ethyl |
| 3. Thionazine | 10. Pirimiphos methyl | 17. Methidathion |
| 4. Ethoprophos | 11. Chlorpyrifos | 18. Ethion |
| 5. Phorate | 12. Malathion | 19. Trithion |
| 6. Fonophos | 13. Parathion | 20. Azinphos methyl |
| 7. Diazinone | 14. Pirimiphos | 21. Azinphos ethyl |

Organophosphorus Pesticides

Column: **MEGA-POF 1 - 0.32mm, 0.25 μ m, 25m**

Catalog Code: C-POF1-032-025-25

Conditions:

Injection: Splitless 240°C.

Detector: FPD 140°C.

Oven Temp.: 100°C (1min), 5.5°C/min, 230°C,
30°C/min, 270°C.

Carrier Gas: Hydrogen, 70kPa.

