


GC Column Dimensions Available

Category	ID	Length	Film Thickness*
Fast	0.05 mm	from 2 to 5 m	0.05 µm 0.10 µm
	0.10 mm	from 2 to 10 m	0.05 µm 0.10 µm 0.20 µm
	0.15 mm	from 2.5 to 30 m	from 0.05 to 0.50 µm
Conventional	0.18 mm 0.20 mm	from 10 to 60 m	from 0.05 µm to 1.20 µm
	0.25 mm	from 10 to 105 m	from 0.05 µm to 3.00 µm
	0.32 mm	from 10 to 105 m	from 0.05 to 5.00 µm
Wide-Bore	0.53 mm	from 10 to 105 m	up to 5.00 µm

Completely customize your product, selecting all combinations of sizes and even asking for out of catalog measures. We can develop ad-hoc solutions for your specific analytical problem without any additional cost.

All our stationary phases are available for Fast, Conventional, Wide-Bore and Multidimensional-GC.



*: the maximum film thickness depends on the stationary phase type